FixNip™ NRI

INSTRUCTIONS FOR USE

Table of Contents:

ENGLISH 2

FRANÇAIS..... 9

ISO SYMBOLS.....17

ENGLISH

FixNip™ NRI - Nipple Reconstruction Implant

PRODUCT DESCRIPTION

The FixNip NRI™ - Nipple Reconstruction Implant, is an implantable hypodermic, silicone implant specially designed for nipple reconstruction following mastectomy. It is comprised of a “floral-shaped” nitinol frame with a silicone capping which allows for: (1) anchoring of the implant to the adjacent fat tissue, (2) providing a softer feel and (3) encapsulating the nitinol frame (Figure 1, Figure 2).

The FixNip NRI™ intended to be implanted within the subcutaneous tissue.

The implant is currently provided in a single size (“one size fits all”).

Figure 1: The FixNip™ – Nipple Reconstruction Implant (NRI)

INTENDED USE

FixNip NRI™ is intended for nipple reconstruction in patients who have undergone mastectomy with removal of the nipple or nipple-areolar complex (NAC).
FixNip NRI™ should be used only by trained professional physicians who are familiar with nipple reconstruction procedures.
FixNip NRI™ is supplied sterile and is intended for one-time use.

CONTRAINDICATIONS

- FixNip NRI is made of Nitinol (NiTi) and contraindicated for patients with known sensitivity to Nickel/Titanium.
 Note: FixNip NRI may contain residuals of Chromium Materials in trace amounts that may cause an allergic reaction
- Tissue thickness of the treated breast in the nipple area is less than 1.0 cm.
- Patient with ongoing infection.
- Pregnant or lactating patient.

WARNINGS 

- The FixNip NRI must be used according to its instructions for use
- The implant is designed for single use and for single patient use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- The implant is sterile. If the package is dry, unopened and undamaged. Do not use if the package is compromised in any way.
- Discard the implant if mishandling has caused possible damage to the device structure or contamination, or if the implant is past its expiration date.
- The device should be used only by a trained physician familiar with the surgical technique for nipple reconstruction.
- The safety and performance of the FixNip NRI when used in conjunction with different breast implant types and/or supportive materials (e.g., mesh) have not been thoroughly evaluated due to limited clinical evidence.
- The implant should be placed in maximum possible contact with healthy, well-vascularized tissue to encourage cellular in-growth and tissue remodeling, at a minimum depth of 0.5 cm.
- The device should not be implanted in infected or potentially infected tissue beds, or over open cavities because infection, migration, or extrusion may result.
- Patients with compromised immune system (such as those with autoimmune disease, diabetes, or undergoing chemotherapy or radiation therapy) may not experience normal wound healing.
- Excessive internal pressure within the implant pocket/reconstructed nipple may increase the risk of extrusion, and local ischemia.

- No studies have been conducted to evaluate the performance of the implant in patients who are pregnant, may become pregnant, or are breast feeding.
- For patients who had been irradiated in the same breast/chest region, consider waiting to implant the device for at least 1 year following their final treatment.
- Tissue at the implant site should be adequately thick (at least 1.0 cm), soft and pliable before implanting the device.
- Patients should wait at least 6 months before getting the treated nipple tattooed. Timing of tattoo placement, quality of the tattooing instruments, and tattooing technique may affect nipple projection.
- Mammography should not be performed on patients treated with the FixNip NRI within 2 months of implantation.
- Suturing the implant in the nipple structure area may result in impaired blood supply, tissue ischemia, necrosis and/or implant exposure. Use only the dedicated holes at the implant base.
- Smoking may interfere with the tissue healing process.
- Lactation may not be feasible after FixNip NRI implantation.

POTENTIAL COMPLICATIONS

Complications that can occur during nipple reconstruction include, but are not limited to:

- Allergic reaction
- Inflammation
- Implant migration
- Implant extrusion/exposure
- Seroma
- Hematoma
- Numbness
- Paresthesia
- Infection
- Tissue ischemia
- Localized necrosis
- Epidermolysis
- Skin discoloration
- Induration
- Implant exposure
- Bleeding
- Insufficient or excessive augmentation
- Pressure sores
- Skin rupture
- Capsular contracture
- Implant inversion (flip over)
- Implant folding
- Discomfort / unnatural feeling
- Pain

INSTRUCTIONS FOR USE

1. Determine local tissue thickness and verify adequacy using a pinch test.
2. Mark exact location of implant placement with round marking ink pen
3. Aseptically prepare the patient and surgical site.
4. Aseptically remove the implant from its pouch and place it in the sterile field. Discard the implant if sterility was compromised.


Figure 2: The FixNip NRI Dimensions

NIPPLE RECONSTRUCTION PROCEDURE

1. Make a 15-25 mm incision a few mm from the intended margin of the reconstructed NAC (Nipple Areola Complex).
2. The FixNip NRI should be placed 3-6 mm from the incision.
3. Develop a 5 mm depth flap in the subcutaneous plane.
4. Perform meticulous hemostasis in the pocket.
5. Fold the FixNip NRI as seen in Figure 3.
6. Insert the FixNip NRI.
 Place the implant. Assure that the implant is flat with no folds in the dissected pocket.
 NOTE: The implant pocket should be developed to allow for exact match in pocket-implant sizes.
 NOTE: To avoid implant migration/displacement it is possible to suture it to the adjacent breast tissue. The suturing is allowed solely through the designated holes located at the outer ring of the implant (see Figure 2).

NOTE: Breast reconstruction by autologous tissue should be well healed before placing the FixNip NRI to ensure the tissue is well vascularized and has good lymphatic drainage

7. After reinforcement, close the incision line in layers in a usual fashion.
 NOTE: It is advisable to dress the reconstructed nipple with soft non-adherent dressing.


Figure 3: Device Implementation Procedure

POSTOPERATIVE CARE

1. Basic local wound care should be applied to the surgical site in the immediate post-operative period.
2. Based on visualized tissue healing, Incision site sutures may be removed at approximately 7-10 days based on visualized tissue healing.
3. Once the skin incision has healed the skin may be tattooed to resemble skin color of a nipple areolar complex. Patients should wait at least 6 months before undergoing tattooing procedure.

STORAGE & EXPIRATION

- The FixNip NRI should be stored in a clean, dry location at room temperature.
- "Use by" date is specified on the implant package.

DISPOSAL

Disposal should be handled in an environmentally sustainable manner according to local regulations. Hazardous waste from contaminated devices or sharps should be disposed of in appropriate containers that meet specific technical requirements.

MRI Safety Information 

The FixNip NRI is MR Conditional. A patient with this Implant may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient:

1.5-Tesla and 3-Tesla MRI

- Static magnetic field of 1.5 or 3 Tesla only.
- Maximum spatial gradient magnetic field of 4000 Gauss/cm (40 T/m) or less
- Maximum MR system reported whole body averaged specific absorption rate (SAR) of <2 W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode
- The presence of the other implants or the health state of the patient may require a reduction in the MR limits.

Under the scan conditions defined above, the FixNip NRI is expected to produce a maximum temperature rise of 1.5°C after 15 minutes of continuous scanning.
In non-clinical testing, the image artifact caused by the device extends approximately 5mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

7-Tesla MRI

- Static magnetic field of 7-Tesla, only
- Maximum spatial gradient magnetic field of 5,000-Gauss/cm (50-T/m)

IMPORTANT NOTE: The 7-Tesla MR system (MAGNETOM, Terra, Siemens Healthineers) approved by the Food and Drug Administration currently has a transmit/receive RF head and a transmit/receive RF knee coil. There is no transmit/receive RF body coil. The FixNip NRI is not implanted in the anatomic regions of the head or the knee.

Therefore, the MRI testing performed on this product was limited to an assessment of magnetic field interactions and artifacts, only. Once a transmit/receive RF body coil becomes commercially available for this application, it will be necessary to conduct further MRI testing to evaluate MRI-related heating for this Implant.

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FixNip™ Ltd. All Rights Reserved 2026
FixNip Ltd.